

TECHNICAL DATA SHEET – Plastic Thinner

PRODUCT: Plastic Thinner “Honey”
TECHNICAL CALLS: 1-800-321-0672



DESCRIPTION:

Plastic Thinner “Honey” is a polyester resin formulated to thin and revitalize body fillers, restoring spreadability and flexibility in just minutes. Eliminates waste.

PART NUMBER:

- 40010 Plastic Thinner “Honey” 16 oz. Can 12 cans / case 14 lbs / case

PRODUCT USES:

Restore older, dried-out body fillers for improved flow-out and spreadability.



MIXING:

Add a small amount of Plastic Thinner to body filler (prior to adding hardener) and continue to add until you reach the desired smoothness and consistency.

APPLICATION:

Apply body filler according to directions on filler’s technical data sheet.



TECHNICAL INFORMATION:

Appearance as Packaged:	Amber Brown
Viscosity @ 77 F:	450 cps (Average)
Density:	9.5 lb/gal (Average)
Catalyst Required:	Benzoyl Peroxide
Catalyzation Ratio:	4.0%
Gel Time @ 77 F:	4.0 – 5.0 minutes
Exotherm Temperature:	220°F (Average)
Exotherm Peak Time:	9.0 minutes (Average)
Shore “D” Hardness Values @ 24 hours:	15-20

ASSOCIATED MSDS: “Plastic Thinner-40010”



HEALTH & SAFETY:

Read all warnings, first aid and safety for all components before using. Keep out of reach of children and animals. Protect hands with impervious rubber gloves. Wear face, skin and eye protection. When sanding, we recommend the use of a respiratory covering device to protect from dust (MSA mask P/N 459029 with MSA cartridge 464029 or equivalent). When using power equipment, refer to power tool manufacturer’s recommendations for safety equipment. USC products are for industrial use by trained professionals only.

Emergency Medical or Spill Control Information:

In U.S. and Canada call CHEMTREC at 1-800-424-9300